INSTRUCTIONS FOR USE

Puritan[®] Sputum Solution

CONTENTS

English.....	2
Français.....	5
Deutsch.....	8
Italiano.....	11
Español.....	14
Svenska.....	17
Nederlandse.....	20

 **Puritan Medical Products Co.**
31 School Street, P.O. Box 149
Guilford, Maine 04443 USA

CALL 800.321.2313 • 207.876.3311
EMAIL sales@puritanmedproducts.com
VISIT puritanmedproducts.com

EC REP EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands 

CH REP MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

Puritan® Sputum Solution

INTENDED USE

Puritan® Sputum Solution is intended for use in the collection and liquefying of clinical sputum specimens for the isolation of bacteria and fungi that cause respiratory tract infections.

SUMMARY & PRINCIPLES

Lower respiratory tract infections are among the most common infectious diseases worldwide.¹ Liquefaction and homogenization of sputum specimens have been recommended to ensure accurate, representative bacterial and fungal cultures.² Puritan Sputum Solution is a nonnutritive phosphate buffered solution containing dithiothreitol (DTT) to act as a liquefying agent for sputum specimens prior to culture.

REAGENTS

Sodium Phosphate, monobasic	Sodium Chloride
Sodium Phosphate, dibasic	DL-Dithiothreitol (DTT)

PRECAUTIONS

- For *in vitro* diagnostic use.
- For single use only.
- For use by trained qualified professionals using aseptic techniques.
- The sterile pipette is to be used for transferring biological specimens (sputum, mucus) from a container into the Sputum Solution. It must not be used for direct collection from the patient.
- Clinical specimens are considered biohazard and must be handled in manner to protect laboratory personnel.
- Clinical samples may contain human pathogens including hepatitis virus and Human Immunodeficiency Virus. Institutional and universally recognized guidelines should be followed when handling items contaminated with blood and other body fluids.³
- Sterilize the unit after use and dispose of it according to biohazard waste disposal regulations.
- Do not use beyond expiry date.
- Do not use if the vial is damaged or there is evidence of contamination, discoloration or leakage.
- Do not ingest the medium.
- The use of this product in combination with diagnostic kits or instruments must be validated by the user prior to use.

STORAGE

For optimum performance, store at 2-25°C (36-77°F). Avoid freezing and excessive heat.^{4,5} The packaging of the Sputum Solution comprises of barrier plastic film together with a modified atmosphere inside the package to maintain stability of the product. The Sputum Solution tubes should be used immediately or stored at 2-8°C and used within 24 hours after opening the package.

MATERIALS PROVIDED

Puritan Sputum Solution is available in product configurations indicated in the table below.

Item No.	Description	Pack Size
SP-100	Natural polypropylene screw-cap tube with 1 mL of Sputum Solution	30 (30/cs)
SP-100 P	Natural polypropylene screw-cap tube with 1 mL of Sputum Solution and one sterile Pasteur pipette for sputum sample transfer	20 (20 kits/cs)

MATERIALS REQUIRED BUT NOT PROVIDED

Materials for the microscopic examination, cultivation, differentiation, and isolation of bacteria from clinical specimens are not provided.

DIRECTIONS FOR USE

Sputum Solution with Sterile Pipette

1. Peel open one sterile pipette.
2. Open the plastic pack containing the vials of Sputum Solution and remove one vial.
3. Insert the pipette into the container of sputum specimen. Squeeze the bulb to obtain approximately 1mL of sputum.
4. Remove the cap of a vial and transfer the contents of the pipette into the Sputum Solution.
5. Replace cap on the Sputum Solution.
6. Leave the inoculated Sputum Solution at room temperature for at least 15 minutes. A prolonged contact time of up to 6 hours does not impair the survival of the bacteria or fungi in the specimen.
7. Vortex the inoculated Sputum Solution for 30 seconds at 2000–2500 rpm or until the mixture has liquefied.
8. Aseptically remove aliquots of the liquefied sputum specimen and inoculate onto appropriate bacteriology culture medium.

SPECIMEN COLLECTION AND HANDLING

Proper specimen collection is critical for successful isolation and identification of infectious organisms. For specific guidance regarding specimen collection procedures, consult published reference manuals.^{6,7}

QUALITY CONTROL

All batches of Puritan Sputum Solution are tested prior to release for pH and further evaluated for their ability to promote growth for the following organisms.

Control	Incubation	Results
<i>Pseudomonas aeruginosa</i> ATCC 27853	Aerobic, 6 h@ Room Temperature	Good Recovery
<i>Haemophilus influenzae</i> ATCC 10211	Aerobic, 6 h@ Room Temperature	Good Recovery
<i>Staphylococcus aureus</i> ATCC 6538	Aerobic, 6 h@ Room Temperature	Good Recovery
<i>Streptococcus pneumoniae</i> ATCC 6305	Aerobic, 6 h@ Room Temperature	Good Recovery
<i>Streptococcus pyogenes</i> ATCC 19615	Aerobic, 6 h@ Room Temperature	Good Recovery
<i>Moraxella catarrhalis</i> ATCC 25238	Aerobic, 6 h@ Room Temperature	Good Recovery
<i>Candida albicans</i> ATCC 10231	Aerobic, 6 h@ Room Temperature	Good Recovery

REPORTING OF SERIOUS INCIDENTS

Any serious incident that has occurred in relation to the device, the patient and/ or the end user should be reported to the manufacturer and, where applicable, the competent authority of the Member State in which the user/ of patient is established.

SYMBOL DEFINITIONS

Consult www.puritanmedproducts.com/symbols-glossary for definition of symbols used in Puritan labeling.

LIMITATIONS

- Puritan Sputum Solution is not tested for acid resistant bacilli culture.
- Performance testing was conducted using laboratory strains spiked into the Sputum Solution tube. Human specimens were not used.

REFERENCES

1. Carroll, K. C. 2002. Laboratory Diagnosis of Lower Respiratory Tract Infections: Controversy and Conundrums. *Journal of Clinical Microbiology*. 40(9): 3115-3120.
2. Hammerschlag, M.R., L. Harding, A. Macone, A. L. Smith, D. A. Goldmann. 1980. Bacteriology of Sputum in Cystic Fibrosis: Evaluation of Dithiothreitol as a Mucolytic Agent. *Journal of Clinical Microbiology*. 11(6): 552-557.
3. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risk related exposure to biological agents at work. *Official Journal of the European Communities*. L 262/21-45.
4. Versalovic, J., K.C. Carroll, G. Funke, J.H. Jorgensen, M.L. Landry, D.W. Warnock. 2011. *Manual of Clinical Microbiology*, 10th ed. American Society for Microbiology. Washington, DC.
5. Miller, J.M. 1996. *A guide to specimen management in clinical microbiology*. American Society for Microbiology. Washington, DC.
6. Forbes, B.A., D.F. Sahm, A.S. Weissfeld. 2007. *Diagnostic Microbiology* 12th ed. Mosby. St. Louis, MO.
7. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, R.H. Tenover. 2003. *Manual of Clinical Microbiology*, 8th ed. American Society for Microbiology. Washington, DC.

Puritan Medical Products Co.
31 School Street, P.O. Box 149
Guilford, Maine 04443 USA

CALL 800.321.2313 • 207.876.3311
EMAIL sales@puritanmedproducts.com
VISIT puritanmedproducts.com

EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

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LITERATUR

1. Carroll, K. C. 2002. Laboratory Diagnosis of Lower Respiratory Tract Infections: Controversy and Conundrums. *Journal of Clinical Microbiology*. 40(9): 3115-3120.
2. Hammerschlag, M.R., L. Harding, A. Macone, A. L. Smith, D. A. Goldman. 1980. Bacteriology of Sputum in Cystic Fibrosis: Evaluation of Dithiothreitol as a Mucolytic Agent. *Journal of Clinical Microbiology*. 11(6): 552-557.
3. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risk related exposure to biological agents at work. *Official Journal of the European Communities*. L 262/21-45.
4. Versalovic, J., K.C. Carroll, G. Funke, J.H. Jorgensen, M.L. Landry, D.W. Warnock. 2011. *Manual of Clinical Microbiology*, 10th ed. American Society for Microbiology. Washington, DC.
5. Miller, J.M. 1996. *A guide to specimen management in clinical microbiology*. American Society for Microbiology. Washington, DC.
6. Forbes, B.A., D.F. Sahm, A.S. Weissfeld. 2007. *Diagnostic Microbiology* 12th ed. Mosby. St. Louis, MO.
7. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, R.H. Tenover. 2003. *Manual of Clinical Microbiology*, 8th ed. American Society for Microbiology. Washington, DC.

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