INSTRUCTIONS FOR USE

Puritan[®] Cary-Blair Medium

CONTENTS

English.....**2**
Français.....**5**
Deutsch.....**8**
Italiano.....**11**
Español.....**14**
Svenska.....**17**
Nederlandse.....**20**

 **Puritan Medical Products Co.**
31 School Street, P.O. Box 149
Guilford, Maine 04443 USA

CALL 800.321.2313 • 207.876.3311
EMAIL sales@puritanmedproducts.com
VISIT puritanmedproducts.com

EC REP EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands 

CH REP MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

Puritan® Cary-Blair Medium

INTENDED USE

Puritan® Cary-Blair is intended for use in the preservation of clinical fecal and rectal specimens containing enteric bacteria from the collection site to the testing laboratory for examination and culture.

SUMMARY & PRINCIPLES

Foodborne illnesses and other diarrheal infections present as a major public health issue. While enteric infections can be caused by different types of bacteria, most routine stool cultures are used to screen for *Salmonella* spp., *Shigella* spp., and *Campylobacter* spp. Cultures for *Vibrio* spp., *Yersinia* spp., and *E. coli* O157:H7 require additional media or incubation conditions and therefore require more extensive preparation.¹

Modified Cary-Blair medium is a nonnutritive balanced salt solution containing phosphates to provide buffering capability, chloride salts to provide essential ions that help maintain osmotic balance. Agar is a solidifying agent and gives a semisolid texture to the medium. Sodium salts and L-cysteine provides a reduced environment. It is recommended for maintaining the viability of enteric bacteria during the transport to the laboratory.²

REAGENTS

Approximate modified Cary-Blair formulation per liter.

Sodium chloride	5.0g	L-cysteine	1.0g
Disodium phosphate	1.1g	Bacteriological agar	1.5g
Sodium thioglycolate	1.5g	Deionized water	1 liter
Calcium chloride	0.09g		

PRECAUTIONS

- For *in vitro* diagnostic use.
- For single use only.
- To be used by trained and qualified professionals using aseptic technique.
- All clinical specimens may contain infectious microorganisms and should be considered biohazards and handled with care. Appropriate personal protective equipment should be worn. Follow laboratory and biosafety guidelines when handling clinical specimens.³⁻⁶
- Refer to the recommendations of the Center for Disease Control and Prevention's Biosafety in Microbiological and Biomedical Laboratories.³⁻⁶
- Sterilize the unit after use and dispose of it according to biohazard waste disposal regulations.
- Do not use beyond expiry date.
- Do not use if the vial is damaged or there is evidence of contamination, discoloration or leakage.
- Do not ingest the medium.
- The use of this product in combination with diagnostic kits or instruments must be validated by the user prior to use.

STORAGE

For optimum performance, store at 2-25°C (36-77°F).

MATERIALS PROVIDED

Puritan Cary-Blair Medium is available in product configurations indicated in the table below.

Item No.	Description	Pack Size
CB-200	Orange polypropylene screw-cap vial with 2 mL Cary-Blair Medium	50 (50/cs)
CB-500	Orange polypropylene screw-cap vial with 5 mL Cary-Blair Medium	50 (50/cs)

MATERIALS REQUIRED BUT NOT PROVIDED

Materials for the microscopic examination, cultivation, differentiation, and isolation of bacteria from clinical specimens are not provided.

DIRECTIONS FOR USE

1. Collect specimen from the rectum or from fresh stool.
2. Remove vial cap and inoculate Cary-Blair medium with rectal or fecal swab specimen.
3. Replace the vial cap, securing tightly. Record patient information in the space provided on the vial label and transport the specimen to the laboratory.

SPECIMEN CULTURES IN THE LABORATORY

In the laboratory, specimens should be processed for bacteriological culture using recommended culture media and laboratory techniques which will depend on the specimen type and the organism under investigation. For recommended culture media and techniques for the isolation and identification of bacteria from clinical swab specimens refer to published microbiology manuals and guidelines.⁷⁻¹¹

QUALITY CONTROL

Each lot of Puritan Cary-Blair medium is tested for sterility, pH, and nonviable bio-burden levels. Representative samples of each lot are further evaluated for their ability to maintain the viability of selected bacterial agents over pre-defined time periods. All bacterial test isolates and testing procedures were established using criteria outlined in the Clinical and Laboratory Standards Institute's M40-A2 document.¹⁰

REPORTING OF SERIOUS INCIDENTS

Any serious incident that has occurred in relation to the device, the patient and/or the end user should be reported to the manufacturer and, where applicable, the competent authority of the Member State in which the user/or patient is established.

SYMBOL DEFINITIONS

Consult www.puritanmedproducts.com/symbols-glossary for definition of symbols used in Puritan labeling.

LIMITATIONS

1. For optimal recovery of *C. difficile*, fecal specimens should be refrigerated at 2-8°C and processed within 48 hours or stored at room temperature (20-25°C) and processed within 24 hours.
2. Reliable specimen collection and transport depends on many factors, including collection and handling, specimen condition, volume, and timing. Best results are achieved when specimens are processed shortly after the time of collection. For detailed information, refer to corresponding reference standards and procedures for optimum collection techniques.^{8, 9, 11-13}
3. Puritan Cary-Blair medium is recommended for the preservation of bacteriological samples only. Viruses, chlamydia, mycoplasma, and ureaplasma require a transport medium formulated specifically for use with these organisms.^{2, 13}
4. Extreme temperature should be avoided during transportation of Puritan Cary-Blair medium.

REFERENCES

1. Centers for Disease Control and Prevention. 2004. Diagnosis and Management of Foodborne Illnesses. Morbid Mortal Weekly Rep. 53: 1-33.
2. Zimbro, M.J., D.A. Power, S.M. Miller, G.E. Wilson, J.A. Johnson. 2009. Difco & BBL Manual of Microbiological Culture Media, 2nd ed. Becton, Dickinson and Company. Sparks, MD.
3. Sewell, D.L. 1995. Laboratory-associated infections and biosafety. Clin. Microbiol. Rev 8:398-405. American Society for Microbiology. Washington, DC.
4. Code of Federal Regulations, title 42, part 72. Interstate shipment of etiologic agents.
5. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risk related exposure to biological agents at work. Official Journal of the European Communities. L 262/21-45.
6. Centers for Disease Control and Prevention. 2009. Biosafety in Microbiological and Biomedical Laboratories, 5th ed. U.S. Department of Health and Human Services, HHS Publication No. (CDC) 21-1112, rev. December 2009.
7. Jorgensen, J.H., M.A. Pfaller, K.C. Carroll, G. Funke, M.L. Landry, S.S. Richter. D.W. Warnock. 2015. Manual of Clinical Microbiology, 11th ed. American Society for Microbiology. Washington, DC.
8. Miller, J.M. 1999. A guide to specimen management in clinical microbiology. American Society for Microbiology. Washington, DC.
9. Forbes, B.A., D.F. Sahm, A.S. Weissfeld. 2007. Diagnostic Microbiology 12th ed. Mosby. St. Louis, MO.
10. CLSI. Quality Control of Microbiological Transport Systems; Approved Standard—Second Edition. CLSI document M40-A2. Wayne, PA: Clinical Laboratory Standards Institute; 2014.
11. Isenberg, H.D. 1998. Collection, Transport and Manipulation of Clinical Specimens. In Essential Procedures for Clinical Microbiology, Ch. 14.12:14-21, 24-27. American Society for Microbiology. Washington, DC
12. Human, R.P., G.A. Jones. 2004. Evaluation of swab transport systems against a published standard. J. Clin. Pathol. 57:762-763 doi:10.1136/jcp.2004.016725.
13. Wasfy, M., B. Oyofu, A. Elgindy, A. Churilla. 1995. Comparison of Preservation Media for Storage of Stool Samples. Journal of Clinical Microbiology. 33(8): 2176-2178.

Puritan[®]
Quality since 1919

Puritan Medical Products Co.
31 School Street, P.O. Box 149
Guilford, Maine 04443 USA

CALL 800.321.2313 • 207.876.3311
EMAIL sales@puritanmedproducts.com
VISIT puritanmedproducts.com

EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

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**BIBLIOGRAFIA**

1. Centers for Disease Control and Prevention. 2004. Diagnosis and Management of Foodborne Illnesses. Morbid Mortal Weekly Rep. 53: 1-33.
2. Zimbro, M.J., D.A. Power, S.M. Miller, G.E. Wilson, J.A. Johnson. 2009. Difco & BBL Manual of Microbiological Culture Media, 2nd ed. Becton, Dickinson and Company. Sparks, MD.
3. Sewell, D.L. 1995. Laboratory-associated infections and biosafety. Clin. Microbiol. Rev 8:398-405. American Society for Microbiology. Washington, DC.
4. Code of Federal Regulations, title 42, part 72. Interstate shipment of etiologic agents.
5. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risk related exposure to biological agents at work. Official Journal of the European Communities. L 262/21-45.
6. Centers for Disease Control and Prevention. 2009. Biosafety in Microbiological and Biomedical Laboratories, 5th ed. U.S. Department of Health and Human Services, HHS Publication No. (CDC) 21-1112, rev. December 2009.
7. Jorgensen, J.H., M.A. Pfaller, K.C. Carroll, G. Funke, M.L. Landry, S.S. Richter. D.W. Warnock. 2015. Manual of Clinical Microbiology, 11th ed. American Society for Microbiology. Washington, DC.
8. Miller, J.M. 1999. A guide to specimen management in clinical microbiology. American Society for Microbiology. Washington, DC.
9. Forbes, B.A., D.F. Sahm, A.S. Weissfeld. 2007. Diagnostic Microbiology 12th ed. Mosby. St. Louis, MO.
10. CLSI. Quality Control of Microbiological Transport Systems; Approved Standard—Second Edition. CLSI document M40-A2. Wayne, PA: Clinical Laboratory Standards Institute; 2014.
11. Isenberg, H.D. 1998. Collection, Transport and Manipulation of Clinical Specimens. In Essential Procedures for Clinical Microbiology, Ch. 14.12:14-21, 24-27. American Society for Microbiology. Washington, DC
12. Human, R.P., G.A. Jones. 2004. Evaluation of swab transport systems against a published standard. J. Clin. Pathol. 57:762-763 doi:10.1136/jcp.2004.016725.
13. Wasfy, M., B. Oyofe, A. Elgindy, A. Churilla. 1995. Comparison of Preservation Media for Storage of Stool Samples. Journal of Clinical Microbiology. 33(8): 2176-2178.

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