INSTRUCTIONS FOR USE

# Puritan<sup>®</sup> Selenite Broth Transport Medium

**CONTENTS**

|                  |           |
|------------------|-----------|
| English.....     | <b>2</b>  |
| Français.....    | <b>5</b>  |
| Deutsch.....     | <b>8</b>  |
| Italiano.....    | <b>11</b> |
| Español.....     | <b>14</b> |
| Svenska.....     | <b>17</b> |
| Nederlandse..... | <b>20</b> |

 **Puritan Medical Products Co.**  
31 School Street, P.O. Box 149  
Guilford, Maine 04443 USA

**CALL** 800.321.2313 • 207.876.3311  
**EMAIL** [sales@puritanmedproducts.com](mailto:sales@puritanmedproducts.com)  
**VISIT** [puritanmedproducts.com](http://puritanmedproducts.com)

**EC REP** EMERGO EUROPE  
Westervoortsedijk 60  
6827 AT Arnhem  
The Netherlands

**CH REP** MedEnvoy Switzerland  
Gotthardstrasse 28  
6302 Zug  
Switzerland

# Puritan® Selenite Broth Transport Medium

## INTENDED USE

Puritan® Selenite Broth Transport Medium is a selective enrichment medium used for the isolation of *Salmonella* spp and *Shigella* spp.

## SUMMARY & PRINCIPLES

*Salmonella* is an important bacterial pathogen of foodborne illness, ranking just behind *C. jejuni* in its frequency.<sup>1</sup> Selenite Broth allows for enhanced growth of *Salmonella* spp in fecal specimens since the pathogen usually represents only a small percentage of the intestinal flora. The peptone provides essential nitrogenous and carbon compounds. Lactose and sodium phosphate maintain a neutral pH. Sodium selenite inhibits many species of gram-positive and gram-negative bacteria including enterococci and coliforms.<sup>2</sup>

## REAGENTS

Formulation per liter.

|                             |                  |
|-----------------------------|------------------|
| Pancreatic Digest of Casein | Sodium Phosphate |
| Lactose                     | Deionized Water  |
| Sodium Selenite             |                  |

pH 7.3 + 0.2 @ 25°C

## PRECAUTIONS

- For *in vitro* diagnostic use.
- For single use only.
- To be used by trained and qualified professionals using aseptic technique.
- Clinical specimens are considered biohazard and must be handled in manner to protect laboratory personnel.
- Clinical samples may contain human pathogens including hepatitis virus and Human Immunodeficiency Virus. Institutional and universally recognized guidelines should be followed when handling items contaminated with blood and other body fluids.<sup>3</sup>
- Sterilize the unit after use and dispose of it according to biohazard waste disposal regulations.
- Do not use beyond expiry date.
- Do not use if the vial is damaged or there is evidence of contamination, discoloration or leakage.
- Do not ingest the medium.
- The use of this product in combination with diagnostic kits or instruments must be validated by the user prior to use.

## STORAGE

For optimum performance, store at 2-25°C. Avoid freezing and overheating.<sup>4,5</sup>

## MATERIALS PROVIDED

Puritan Selenite Broth Transport Medium is available in product configurations indicated in the table below.

| Item No. | Description   | Pack Size  |
|----------|---|------------|
| SB-500   | White polypropylene screw-cap tube with 5 mL of Selenite Broth Medium | 50 (50/cs) |

## MATERIALS REQUIRED BUT NOT PROVIDED

Materials for the microscopic examination, cultivation, differentiation, and isolation of bacteria from clinical specimens are not provided.

## DIRECTIONS FOR USE

### Selenite Broth Collected Sample

1. Vortex the inoculated Selenite Broth transport medium for approximately 10 seconds.
2. Incubate inoculated Selenite Broth transport medium at 35 ± 2°C for 18-24 hours.
3. After incubation, streak the sample on a surface of specific agar plate by using a swab or removing aliquots of the Selenite Broth transport medium and inoculate on to the agar plate.

### Fecal Opti-Swab® Collected Sample

1. Obtain tubes of Selenite Broth transport medium and unscrew cap.
2. Vortex the inoculated Fecal Opti-Swab® for approximately 10 seconds.
3. Unscrew the cap and aseptically transfer the swab from the Fecal Opti-Swab® to the Selenite Broth transport medium using sterile forceps.
4. Replace cap on both Fecal Opti-Swab® and Selenite Broth.
5. Follow the procedures stated above for Selenite Broth Collected Sample.

## SPECIMEN COLLECTION AND HANDLING

Specimens suitable for culture may be handled using various techniques. For detailed guidance, refer to appropriate references.<sup>6,7</sup> Specimens should be obtained before antimicrobial agents have been administered.

## QUALITY CONTROL

All batches of Puritan Selenite Broth Transport Medium are tested prior to release for pH and further evaluated for their ability to promote growth of *Salmonella* spp and *Shigella* spp, and suppress enterococci and coliforms over predefined time periods. All bacterial test isolates and testing procedures were established using criteria outlined in the Clinical and Laboratory Standards Institute's M22-A3 document and dehydrated media manufacturer recommendations where applicable.<sup>2, 8</sup>

| Control   | Incubation                  | Results    |
|---|-----------------------------|------------|
| <i>Salmonella enterica</i> Serovar Typhimurium ATCC 14028 | Aerobic, 18-24 hr @ 35-37°C | Growth     |
| <i>E. coli</i> ATCC 25922                                 | Aerobic, 18-24 hr @ 35-37°C | Inhibition |
| <i>Shigella sonnei</i> ATCC 9290                          | Aerobic, 18-24 hr @ 35-37°C | Growth     |

## REPORTING OF SERIOUS INCIDENTS

Any serious incident that has occurred in relation to the device, the patient and/ or the end user should be reported to the manufacturer and, where applicable, the competent authority of the Member State in which the user/ of patient is established.

## SYMBOL DEFINITIONS

Consult [www.puritanmedproducts.com/symbols-glossary](http://www.puritanmedproducts.com/symbols-glossary) for definition of symbols used in Puritan labeling.

## LIMITATIONS

Definitive identification of *Salmonella* spp and *Shigella* spp requires additional and/or serological tests. Refer to appropriate reference standards for further instructions.<sup>6,7</sup>

REFERENCES

1. Centers for Disease Control and Prevention. 2004. Diagnosis and Management of Foodborne Illnesses. *Morbidity and Mortality Weekly Report*. 53: 1-33.
2. Zimbro M.J, D.A. Power. 2003. *Difco & BBL Manual: Manual of Microbiological Culture Media*. Becton, Dickinson, and Company. Sparks, MD.
3. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risk related exposure to biological agents at work. *Official Journal of the European Communities*. L 262/21-45.
4. Versalovic, J., K.C. Carroll, G. Funke, J.H. Jorgensen, M.L. Landry, D.W. Warnock. 2011. *Manual of Clinical Microbiology*, 10th ed. American Society for Microbiology. Washington, DC.
5. Miller, J.M. 1996. *A guide to specimen management in clinical microbiology*. American Society for Microbiology. Washington, DC.
6. Forbes, B.A., D.F. Sahm, A.S. Weissfeld. 2007. *Diagnostic Microbiology* 12th ed. Mosby. St. Louis, MO.
7. Murray, P.R., E.G. Baron, J.H. Jorgensen, M.A. Pfaller, R.H. Tenover, R. Tenover. 2003. *Manual of Clinical Microbiology*, 8th. American Society for Microbiology, Washington, DC.
8. CLSI. *Quality Control for Commercially Prepared Microbiological Culture Media; Approved Standard-Third Edition*. CLSI document M22-A3. Wayne, PA. Clinical and Laboratory Standards Institute; 2004.

**Puritan Medical Products Co.**  
31 School Street, P.O. Box 149  
Guilford, Maine 04443 USA

**CALL** 800.321.2313 • 207.876.3311  
**EMAIL** [sales@puritanmedproducts.com](mailto:sales@puritanmedproducts.com)  
**VISIT** [puritanmedproducts.com](http://puritanmedproducts.com)

EMERGO EUROPE  
Westervoortsedijk 60  
6827 AT Arnhem  
The Netherlands

MedEnvoy Switzerland  
Gotthardstrasse 28  
6302 Zug  
Switzerland

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## REFERENSER

1. Centers for Disease Control and Prevention. 2004. Diagnosis and Management of Foodborne Illnesses. Morbid Mortal Weekly Rep. 53: 1-33.
2. Zimbro M.J, D.A. Power. 2003. Difco & BBL Manual: Manual of Microbiological Culture Media. Becton, Dickinson, and Company. Sparks, MD.
3. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risk related exposure to biological agents at work. Official Journal of the European Communities. L 262/21-45.
4. Versalovic, J., K.C. Carroll, G. Funke, J.H. Jorgensen, M.L. Landry, D.W. Warnock. 2011. Manual of Clinical Microbiology, 10th ed. American Society for Microbiology. Washington, DC.
5. Miller, J.M. 1996. A guide to specimen management in clinical microbiology. American Society for Microbiology. Washington, DC.
6. Forbes, B.A., D.F. Sahm, A.S. Weissfeld. 2007. Diagnostic Microbiology 12th ed. Mosby. St. Louis, MO.
7. Murray, P.R., E.G. Baron, J.H. Jorgensen, M.A. Pfaller, R.H. Yolden. 2003. Manual of Clinical Microbiology, 8th. American Society for Microbiology, Washington, DC.
8. CLSI. *Quality Control for Commercially Prepared Microbiological Culture Media; Approved Standard-Third Edition*. CLSI document M22-A3. Wayne, PA. Clinical and Laboratory Standards Institute; 2004.

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EMERGO EUROPE  
Westervoortsedijk 60  
6827 AT Arnhem  
Nederländerna

MedEnvoy Switzerland  
Gotthardstrasse 28  
6302 Zug  
Schweiz

# Puritan® Selenite Broth-transportmedium

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