INSTRUCTIONS FOR USE

Puritan® MRSA Transport Medium

CONTENTS

English.....**2**
Français.....**5**
Deutsch.....**8**
Italiano.....**11**
Español.....**14**
Svenska.....**17**
Nederlandse.....**20**

 **Puritan Medical Products Co.**
31 School Street, P.O. Box 149
Guilford, Maine 04443 USA

CALL 800.321.2313 • 207.876.3311
EMAIL sales@puritanmedproducts.com
VISIT puritanmedproducts.com

EC REP EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands 

CH REP MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

Puritan® MRSA Transport Medium

INTENDED USE

The Puritan® MRSA Transport Medium is an enrichment medium used for isolation of *Staphylococcus aureus* spp from clinical specimens, particularly methicillin-resistant *Staphylococcus aureus* (MRSA).

SUMMARY & PRINCIPLES

S. aureus is one of the most common causes of skin and soft tissue infection in both the health care and community settings.¹ Enrichment broths are commonly used to increase sensitivity testing for MRSA by increasing isolation rates. Tryptic Soy Broth (TSB) serves as the base medium for Puritan MRSA Transport Medium to enhance the growth of *S. aureus*. TSB contains enzymatic digest of casein and enzymatic digest of soybean meal, which provide amino acids and complex nitrogenous compounds that promote microbial growth. Dextrose acts as a carbon energy source that facilitates growth. Dipotassium phosphate acts as a buffering agent. Sodium chloride is added to inhibit or partially inhibit microorganisms other than *S. aureus*.

REAGENTS

Formulation per liter.

TSB powder	30.0g
Sodium Chloride	25.0g
Deionized Water	1000 mL

pH 7.3 + 0.2 @ 25°C

PRECAUTIONS

- For *in vitro* diagnostic use.
- For single use only.
- To be used by trained and qualified professionals using aseptic technique.
- Clinical specimens are considered a biohazard and must be handled in a manner to protect laboratory personnel.
- Clinical samples may contain human pathogens including hepatitis virus and Human Immunodeficiency Virus. Institutional and universally recognized guidelines should be followed when handling items contaminated with blood and other body fluids.²
- Sterilize the unit after use and dispose of it according to biohazard waste disposal regulations.
- Do not use beyond expiry date.
- Do not use if the vial is damaged or there is evidence of contamination, discoloration or leakage.
- Do not ingest the medium.
- The use of this product in combination with diagnostic kits or instruments must be validated by the user prior to use.

STORAGE

For optimum performance, store at 2-25°C. Avoid freezing and overheating.^{3,4}

MATERIALS PROVIDED

Puritan MRSA Enrichment Broth is available in product configurations indicated in the table below.

Item No.	Description	Pack Size
MRSA-200	Pink polypropylene screw-cap vial with 2 mL MRSA Medium	50 (50/cs)

MATERIALS REQUIRED BUT NOT PROVIDED

Materials for the microscopic examination, cultivation, differentiation, and isolation of bacteria from clinical specimens are not provided.

DIRECTIONS FOR USE

MRSA Transport Medium Collected Sample

1. Vortex the inoculated MRSA transport medium for approximately 10 seconds.
2. Incubate inoculated MRSA transport medium at 35 ± 2°C.
3. Examine the MRSA transport medium for growth after 18-24 hours.
4. Aseptically remove aliquots of the MRSA transport medium and inoculate on to an appropriate selective agar plate.

Opti-Swab® Liquid Amies Collected Sample

1. Obtain tubes of MRSA Transport Medium and unscrew cap.
2. Vortex the inoculated Opti-Swab® Liquid Amies for approximately 10 seconds.
3. Unscrew the cap and aseptically transfer the swab from the Opti-Swab® Liquid Amies to the MRSA Transport Medium using sterile forceps.
4. Replace cap on both Opti-Swab® Liquid Amies and MRSA Transport Medium.
5. Follow the procedures stated above for MRSA Transport Medium Collected Sample.

SPECIMEN COLLECTION AND HANDLING

Specimens suitable for culture may be handled using various techniques. For detailed guidance, refer to appropriate references.^{5,6} Specimens should be obtained before antimicrobial agents have been administered.

QUALITY CONTROL

All batches of Puritan MRSA Transport Medium are tested prior to release for pH and further evaluated for their ability to promote growth of the following organisms.

Control	Incubation	Results
Methicillin-resistant <i>Staphylococcus aureus</i> ATCC 43300	Aerobic, 48 hr @ room temperature	Good recovery
<i>Staphylococcus aureus</i> ATCC 6538	Aerobic, 48 hr @ room temperature	Good recovery

REPORTING OF SERIOUS INCIDENTS

Any serious incident that has occurred in relation to the device, the patient and/or the end user should be reported to the manufacturer and, where applicable, the competent authority of the Member State in which the user/ of patient is established.

SYMBOL DEFINITIONS

Consult www.puritanmedproducts.com/symbols-glossary for definition of symbols used in Puritan labeling.

LIMITATIONS

Definitive identification of MRSA requires additional and/or serological tests. Refer to appropriate reference standards for further instructions.^{5,6}

REFERENCES

1. National Health and Nutrition Examination Survey (NHANES). 2000. Specimen Collection Procedures Manual.
2. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risk related exposure to biological agents at work. Official Journal of the European Communities. L 262/21-45.
3. Versalovic, J., K.C. Carroll, G. Funke, J.H. Jorgensen, M.L. Landry, D.W. Warnock. 2011. Manual of Clinical Microbiology, 10th ed. American Society for Microbiology. Washington, DC.
4. Miller, J.M. 1996. A guide to specimen management in clinical microbiology. American Society for Microbiology. Washington, DC.
5. Forbes, B.A., D.F. Sahm, A.S. Weissfeld. 2007. Diagnostic Microbiology 12th ed. Mosby. St. Louis, MO.
6. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, R.H. Tenover. 2003. Manual of Clinical Microbiology, 8th ed. American Society for Microbiology. Washington, DC.

Puritan Medical Products Co.
31 School Street, P.O. Box 149
Guilford, Maine 04443 USA

CALL 800.321.2313 • 207.876.3311
EMAIL sales@puritanmedproducts.com
VISIT puritanmedproducts.com

EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

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