



Quality Medical Products Since 1919

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

An affiliate of Hardwood Products Company LLC

Tel 800-321-2313 (US and Canada)
207-876-3311

Fax 800-323-4153 (US and Canada)
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

CALGISWAB® CALCIUM ALGINATE-TIPPED APPLICATORS

These speciality applicators are tipped with calcium alginate fiber, an organic material derived from seaweed and carded in Puritan's facility to assure conformance with medical device specifications. Calgiswab applicators are available in a variety of tip sizes and shaft materials and dimensions. Commonly used for specimen collection in pediatrics, E.N.T., asymptomatic GC, public health microbiology, ophthalmology, food pathology, and other areas where fine tips, flexible shafts, or biodegradable fibers are desirable.

	Item No.	Description	Shaft Length & Dia.	Packing	Weight
	25-800 A 50	Nasopharyngeal, aluminum shaft, sterile, Type 1	5-1/2" x .035"	10/50/1 per pkg.	4 lbs.
	25-801 A 50	Urethro-genital, aluminum shaft, sterile, Type 4	5-1/2" x .025"	10/50/1 per pkg.	4 lbs.
	25-806 1PA 25-806 2PA	Plastic shaft, sterile, Type 3	6" x 1/10"	10/100/1 per pkg. 10/100/2 per pkg.	6 lbs. 7 lbs.
	25-806 1WA 25-806 2WA	Wood shaft, sterile, Type 2	6" x 1/10"	10/100/1 per pkg. 10/100/2 per pkg.	6 lbs. 7 lbs.
		806-PA 5M	Plastic shaft, Type 3	6" x 1/10"	5/10/100
806-WA 5M		Wood shaft, Type 2	6" x 1/10"	5/10/100	9 lbs.

Products are not shown actual size. Rayon, cotton, or polyester available with aluminum or stainless steel shaft upon request.

Did you know?

Calgiswabs are well-suited for surface testing in the pharmaceutical field. The calcium alginate fiber used in urethra-genital swabs contains no fatty acids that may be inhibitory to many strains of gonorrhea. In addition, the alginate fibers are slowly dissolved by body fluids and tissues, and any fibers that remain after testing will quickly be eliminated.

