

# Cap-Shure®

by Puritan®

**Reorder No. 25-807 1PC EC**

**Pack: 10/50/1**

Sterile

7" Plastic Shaft

Cotton Tip with Clear Aerated Cap

Individually Wrapped

**Reorder No. 25-806 1WC EC**

**Pack: 10/50/1**

Sterile

6" Wood Shaft

Cotton Tip with Clear Aerated Cap

Individually Wrapped

## Better sampling — better bottom line



## More convenient. Enhanced reliability.

Now there's a tip-protected applicator that fits both your sampling needs **and** your budget. It's the Puritan Cap-Shure, designed for convenient, reliable specimen, DNA or evidence collection. Sterile cotton-tipped applicator comes as a unit with a reclosable, ventilated cap for air drying and transporting samples.

Quality performance. Budget-conscious price. From the most trusted name in single-use medical products. For more information **contact your Puritan distributor, or call 1-800-321-2313.** Other tip materials are available upon request. Please call to inquire about custom orders.



Quality Medical Products Since 1919

Puritan Medical Products Company LLC  
Guilford, Maine 04443-0149 U.S.A.

Tel 800-321-2313 (US and Canada)  
207-876-3311

Fax 800-323-4153 (US and Canada)  
207-876-3130

Web [www.puritanmedproducts.com](http://www.puritanmedproducts.com)

Email [sales@puritanmedproducts.com](mailto:sales@puritanmedproducts.com)

ISO 9001:2000 CE

