



# Calgiswab® Calcium Alginate Tipped Applicator

Regular tip, wood handle



- Calcium alginate tip

- Rigid wood handle

Not actual size

## Description

This specialty applicator is tipped with calcium alginate fiber, an organic material derived from seaweed and carded in Puritan's facility to assure conformance with medical device specifications. It is bonded securely to a rigid wooden handle by an aqueous based adhesive.

## Dimensions

**6" Overall Length (152.4 mm)**

### Tip

- 0.219" Dia. (5.563 mm)
- 0.641" Length (16.281 mm)

### Handle

- 0.1" Dia. (2.54 mm)
- 5.812" Length (147.625 mm)

## Qualities

- Medical grade quality
- Calcium alginate tip contains no fatty acids, which may be inhibitory to many strains of bacteria
- Made in USA

## Applications

- Specimen collection for testing and screening
- Wound care

Ref #	Packaging
806-WA 5M Non-sterile	5/10/100 (5 boxes per case, 10 bags per box, 100 swabs per bag)
25-806 1WA Sterile	10/100/1 (10 boxes per case, 100 individually wrapped packages per box, 1 p/pkg)



**Puritan®**