



# Histobrush® Cytology Collection Device

Polystyrene handle



-Soft, tapered nylon bristles

- Polystyrene handle

\*Image not actual size

## Description

The Puritan Histobrush® is made with soft, tapered nylon bristles that is gentle on patients.

## Dimensions

**7" Overall Length (177.8 mm)**

### Tip

- 0.2" Dia. (4.775 mm) at tapered end
- 1" Length (25.4 mm)

### Handle

- 0.121" Dia. (3.073 mm)
- 6.219" Length (157.963 mm)

## Qualities

- Polystyrene handle
- Soft, tapered nylon bristles

## Applications

- Cervical specimen collection
- Buccal cell collection
- DNA collection
- Veterinary use

Ref #	Packaging
2199 Non-sterile	10/100 (10 boxes per case, 100 per box)
25-2199 Sterile	10/100/1 (10 boxes per case, 100 individually wrapped packages per box, 1 p/pkg)



**Puritan®**