



PurFlock Ultra® Flocked Swab

Mini-Tip with Polystyrene Handle - no breakpoint



Description

An ideal choice in flock swabs with excellent absorption and elution of specimens when HydraFlock's high levels of absorbency aren't necessary. PurFlock Ultra's high purity characteristics makes it an ideal swab for various rapid diagnostic tests.

Patented: US, Canada, Europe, Australia, Japan, China
Patent pending: Hong Kong

Dimensions

6.04" Overall Length (153.416 mm)

Tip

- 0.135" Dia. (3.429 mm)
- 0.800" Length (20.32 mm)

Handle

- 0.098" Dia. (2.489 mm) tapered to 0.049" (1.245 mm)
- 6" Length (152.654 mm)

Qualities

- Medical grade quality
- Highly absorbent tip
- Soft smooth tip
- Made in USA

Applications

- For nasopharyngeal and urethral cell collection
- Rapid diagnostic testing
- Specimen collection
- For collection of bacteria or micro-organisms
- Forensic evidence collection

Ref #	Packaging
3316-U Non-Sterile	10/100 (10 bags of 100 swabs per bag)
25-3316-U Sterile	10/50/1 (10 boxes of 50 individually wrapped swabs)

