



PurFlock Ultra® Flocked Swab - Dry Transport System

Mini-tip, polystyrene handle with breakpoint at 80 mm



Description

Puritan's polypropylene dry transport tube and plug-style cap is perfect for specimen collection and transportation to the lab or use in diagnostic kit. PurFlock Ultra® flocked swabs provide excellent absorption and elution of specimens when HydraFlock's high levels of absorbency aren't necessary. PurFlock Ultra's high purity characteristics makes it an ideal swab for various rapid diagnostic tests.

Dimensions

Overall Length

- Tube: 6.062" (153.975 mm)
- Swab: 6.04" (153.416 mm)

Tip

- 0.135" Dia. (3.429 mm)
- 0.800" Length (20.32 mm)

Handle

- 0.098" Dia. (2.489 mm) tapered to 0.049" (1.245 mm)
- 6.0" Length (152.4 mm)

Qualities

- Handle with breakpoint at 80 mm
- Medical grade quality
- Highly absorbent capacity
- Packaged in a reclosable dry transport tube for ease in transport to the lab
- Made in USA

Applications

- For nasopharyngeal and urethral cell collection
- Collection of bacteria or microorganisms
- Forensic evidence collection
- Specimen collection

Ref #	Packaging
25-3316-U EMB 80 BT Sterile	5/100/1 (5 boxes, 100 tubes per box, 1 swab per tube)



Puritan®