




# PurFlock Ultra® Flocked Swab

Elongated Tip with Polystyrene Handle - 100mm breakpoint



## Description

PurFlock's purity make it ideal for PCR, molecular assays, DNA testing, forensic applications as well as direct antigen testing. Its uniquely designed fiber structure accounts for better cell yield in the specimen where sample size is limited.

Patented: US, Canada, Europe, Australia, Japan, China  
Patent pending: Hong Kong

## Dimensions

6" Overall Length (152.4 mm)

### Tip

- 0.219" Dia. (5.56 mm)
- 0.685" Length (17.399 mm)

### Handle

- 0.1" Dia. (2.54 mm) tapered to 0.067 (1.702 mm)
- 5.96" Length (151.384 mm)

## Qualities

- Highly absorbent tip
- Soft smooth tip
- Made in USA

## Applications

- Buccal cell / DNA collection
- Forensic evidence collection
- Rapid diagnostic testing
- Specimen collection

Ref #	Packaging
3706-U Non-Sterile	10/100 (10 bags per case, 100 per bag)
25-3706-U Sterile	10/50/1 (10 boxes per case, 50 individually wrapped swabs, 1 swab p/pkg)



**PurFlock  
ULTRA®**