




Cotton Tipped Applicator

Standard Tip, Wood Handle



Description

This standard cotton-tipped applicator is produced with high quality USP grade cotton. It is bonded securely to a wood handle by an aqueous based adhesive. It is lint and residue free, making it an ideal applicator for cleaning sensitive components. Imported swab.

Dimensions

6" Overall Length (152.4 mm)

Tip

- 0.19" Dia. (4.7 mm)
- 0.60" Length (15 mm)

Handle

- 0.093" Dia. (2.35 mm)
- 5.91" Length (150 mm)

Qualities

- Standard cotton tip
- Sturdy wood handle
- Highly absorbent, does not swell in solvents
- Lintless and residue free
- Low in NVR

Applications

- Cleaning sensitive components
- Cleaning disk drives, optics, computer parts, etc.
- Precision application of solutions or adhesives
- Easy disposal

Ref #	Packaging
896-WC Non-sterile	10/100 (10 bags of 100 swabs per bag = 1,000/ cs)



Puritan®



Cotton Tipped Applicator

Standard Tip, Wood Handle

Ref #: 896-WC

Report of Analysis: NVR/Ion/Silicone Oil

Parameter	Result	Unit	Specification
1. Chloride	0.772	µg/ tip	< 8.5
2. NVR (DI)	0.473	mg/ tip	< or = 0.53
3. NVR (IPA)	0.092	mg/ tip	< or = 0.15
4. NVR (Acetone)	0.288	mg/ tip	< or = 0.55
5. Silicone Oil	-	ng/ tip	< 1.3



Puritan[®]